


1 inch



- Full Face Delivery
- Direct Nasal Delivery
- Nasal Delivery
- Oral Delivery




Infinity™481
DIRECT NASAL MASK



FlexiFit™407
PREMIUM NASAL MASK



FlexiFit™406
PETITE NASAL MASK



FlexiFit™431
FULL FACE MASK

185044349REVA

SleepStyle™
SLEEP PERFORMANCE FOR
AN ENERGIZED LIFESTYLE

Fisher & Paykel
HEALTHCARE

Fisher & Paykel Multi-Mask Sizing Gauge

Instructions

1. Print Document
(Note) do not choose to scale document to paper size when printing.
2. Verify correct scale, by measuring 1 inch marker.
3. Cut out paper inside the mask outline.
4. Choose correct size according to fit (406 & 407) or size lines (431, 432 & 481)