


STANDARD




ALIGN WITH
BOTTOM OF NOSE

LARGE



ALIGN WITH
BOTTOM OF NOSE

SHALLOW




ALIGN WITH
BOTTOM OF NOSE

Follow these two steps to choose the correct size:



1 Check the nostril width for clearance



2 Check the nasal bridge width for clearance.

RESMED

Customer Service 1-800-424-0737 (US only)
www.resmed.com

Mirage and ACTiva are trademarks of ResMed Ltd. © 2003 ResMed Ltd.

MIRAGE
ACTiva
ACTIVECELL TECHNOLOGY

FITTING TEMPLATE