

ResMed Mirage Vista Sizing Gauge

60940/4 04 04

MIRAGE VISTA™
Nasal Mask

FITTING TEMPLATE

RESMED

Customer Service 1-800-424-0737 (US only)
www.resmed.com
Mirage and Vista are trademarks of ResMed Ltd. ©2004 ResMed Ltd.

DEEP
CUT OUT
ALIGN WITH BOTTOM OF NOSE

STANDARD
CUT OUT
ALIGN WITH BOTTOM OF NOSE

SMALL
CUT OUT
ALIGN WITH BOTTOM OF NOSE

1 inch

Follow these two steps to choose the correct size:

1 Check the nostril width for clearance.

2 Check the nasal bridge width for clearance.

INSTRUCTIONS:

1. Print out
2. Measure inch marker to verify scale
3. Cut out grey areas
4. Choose size based on instruction printed on gauge